
## Product portfolio overview



\*For information on compatibility, please visit www.bernafon.com/products/accessories.

Please note that availability of some features may be dependent on technical constraints of hearing aid style. Please check Product Information of desired style for more details.

## User benefits



#### DECS™

Dynamic Environment Control System<sup>™</sup>, or DECS<sup>™</sup>, uses a combination of systems to ensure amplification of sound happens in real time, enhancing speech understanding in challenging listening situations, and providing users with the most realistic sound experience possible.



### Sound

Strong signal processing provides users with better speech understanding in dynamic listening situations and a natural sound experience across all frequencies.



### Speech

The right amount of directionality is applied to offer an optimized binaural signal-to-noise ratio across the various frequencies. By targeting and reducing noise where appropriate, users can concentrate on what's important to them.



### Hearing comfort

Different noise and feedback management systems reduce disturbing noise without adversely affecting speech providing users with more hearing comfort.



### **Operating comfort**

Various accessories and connectivity options enable direct sound streaming into hearing aids for more control and operating comfort.



### Individualization

Various individualization options are available for a more personalized fit.



## Model overview

# non l'

Technical overview	miniRITE	miniRITE T	miniRITE T R	BTE 105	lic	CIC	ITC <sup>1)</sup>	ITE HS <sup>1)</sup>	ITE FS <sup>1)</sup>
Viron 9 7 5	•	•	•	•					
Zerena 9 7 5 3 1	•	•		•	•	•	•	•	•
Battery size	312	312	312+	13	10	10	312	312, 13	312, 13
Fitting level	60-, 85-, 100-, 105-Speaker	60-, 85-, 100-, 105-Speaker	60-, 85-, 100-, 105-Speaker	105-Speaker	75-, 85-Speaker	75-, 85-Speaker	75-, 85-, 90-, 100-Speaker	75-, 85-, 90-, 100-Speaker	75-, 85-, 90-, 100-Speaker
miniFit thin tube				•					
Earhook				•					
Li-ion rechargeable battery			•						
ZPower rechargeable battery	•								
Dual-radio 2.4 GHz technology	•	•	•	•			0	0	0
NFMI technology	•	•	•	•		0	•	•	•
Single push button	•					0	0	0	0
Double push button		•	•	•					
Volume wheel							0	0	0
Telecoil		•	•	•			0	0	0
Auto Telephone (detection) <sup>2)</sup>	•	•		•			0	0	0
Hydrophobic coating	•	•	•	•	•	•	•	•	•
IP rating	IP68	IP68	IP68	IP68	IP68	IP68	IP68	IP68	IP68

Connectivity overview <sup>3)</sup>													
•	•	•	•	•	•	•							
			•										
•	•	•	•	•	•	•							
•	•	•	•	•	•	•							
•	•	•	•	•	•	•							
•	•	•	•	•	•	•							
	• • • • • • • • • • • • • • • • • • • •	• • • • • • • • • • • •											

Included Optional

<sup>1)</sup> ITE hearing aids with battery size 312 can accommodate either 2.4 GHz wireless or telecoil. Choose battery size 13 if both options are required.

<sup>2)</sup> Only available in Zerena hearing aids.

<sup>3)</sup> For hearing aids with dual-radio 2.4 GHz technology.

## BTE color options



## ITE color options



\* Not available for Zerena 3|1

\* For IIC only

Apple, the Apple logo, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.


### Accessories



### SoundClip-A



### RC-A (remote control)



### TV-A (TV adapter)



### EasyControl-A app

